

POCT TESTS

POCT TESTS	PLATFORM	SPECIMEN	TUBE TYPE
Activated Clotting Time	I-STAT	Whole blood arterial or venipuncture	Syringe or tube with NO anticoagulant, clot activator or serum separator
	HEPCON	Whole blood arterial or venipuncture	Syringe with NO heparin or anticoagulant
	HEMOCHRON RESPONSE	Whole blood arterial or venipuncture	Syringe with NO heparin or anticoagulant
Blood Gas / Electrolyte / Hct / Lactate / Glucose	I-STAT	Whole blood arterial or venipuncture	Lithium or sodium heparinized syringe
CHEM 8+	I-STAT	Whole blood venipuncture	Lithium or sodium heparinized syringe
CKMB	I-STAT	Whole blood venipuncture	Lithium or sodium heparinized syringe (must be at least half-filled); if in non-heparinized syringe or tube- sample must be tested within 1 minute
Coagulation – PT/INR	I-STAT	Whole blood capillary and venipuncture	Syringe without anticoagulant or fingerstick
Glucose	NOVA Glucose METER	Whole blood capillary or venipuncture	Fingerstick or Syringe without anticoagulant
Hemoglobin	HEMOCUE	Whole blood capillary, venous or arterial venipuncture	Fingerstick or syringe without anticoagulant
Hemoglobin A1C	DCA2000	Whole blood capillary, venous or arterial venipuncture	Fingerstick or syringe without anticoagulant
Occult Blood	HEMOCCULT	Stool	Slides
Oxygen Saturation	AVOX1000E	Whole Blood or plasma	Lithium or sodium heparinized syringe or vacutainer tubes
Fern Test	MICROSCOPE	Vaginal fluid	Swab
pH	pH TEST STRIP	Nasal or gastric fluid	Swab
KOH Prep	MICROSCOPE	Vaginal fluid	Swab
Wet Prep	MICROSCOPE	Vaginal fluid	Swab
Rapid HIV-1 Antibody Test	ORAQUICK Kit	Oral fluid	Swab
Rapid Strep-A	ACCEAVA KIT	Throat fluid	Swab
Urinary Creatinine/Microalbumin	DCA2000	Fresh urine	Yellow top tube (no additives)
Urinalysis	CLINITEK	Fresh urine	Yellow top tube (no additives)
Urine Pregnancy	ICON25 kit	Fresh urine	Yellow top tube (no additives)
Whole Blood Creatinine	I-STAT	Whole blood venipuncture	Lithium or sodium heparinized syringe
Clotting Factors/Time	TEG	Native whole blood	Syringe