

Icahn School of Medicine at Mount Sinai Autoclave use log

Building	
Floor/Room	
DSO	

Make	
Model	
Serial Number	

Cycle number	User	Date	Material	Process type ¹	Cycle time (minutes)	Heat-sensitive tape ²	Chemical integrator ³	Biological indicator ⁴

¹ Process type: Liquid cycle or Dry cycle
² Results indicated (+) positive or (-) negative for color change
³ Results indicated as "Accept" or "Reject"
⁴ Results indicated (+) positive or (-) negative for color change/turbidity indicators of microbial growth